

## 特殊应用类真空采血管

Special Vacuum Blood Collection Tube



### Labtub® 血液DNA管

Labtub® Blood DNA Tube

用于采集、抗凝、储存、运输血液标本，稳定基因组DNA

Used for blood collection, anticoagulation, storage, transport, stabilization of genomic DNA

- 稳定期限

Stability Period

15-25°C稳定14天，2-8°C稳定28天，35°C时至少稳定3天，-20°C至少稳定6个月

14 days at 15-25°C, 28 days at 2-8°C, at least 3 days at 35°C, at least 6 months at -20°C

- 标识

Identification

浅紫色胶塞/透明安全头盖，标签印有不同编号的二维码

Light purple rubber stopper / transparent safety cap, unique number of 2-D barcode on label

产品代码 No.	材质 Material	采血管尺寸 Size	公称容量 Nominal Capacity
25022431	PET	13*100mm	5ml
20023431	PET	16*100mm	10ml



## Labtub® 血液RNA管 Labtub® Blood RNA Tube

用于采集、存储、运输血液标本，稳定基因组RNA  
Used for blood collection, anticoagulation, storage, transport, stabilization of genomic RNA

- 稳定期限**  
**Stability Period**  
 15-25°C至少3天，2-8°C至少5天，-20- -70°C至少8年  
 At least 3 days at 15-25°C, at least 5 days at 2-8°C, at least 8 years at -20- -70°C
- 标识**  
**Identification**  
 浅蓝色胶塞/透明安全头盖，标签印有不同编号的二维码  
 Light blue rubber stopper / transparent safety cap, unique number of 2-D barcode on label

产品代码 No.	材质 Material	采血管尺寸 Size	添加剂体积 Additive Volume	公称容量 Nominal Capacity
29023441	PET	16*100mm	6ml	3ml



## Labtub® 血液ccfDNA管 Labtub® Blood ccfDNA Tube

Stabilization of circulating, cell-free DNA ( ccfDNA )

用于采集、抗凝、储存、运输血液标本，稳定循环游离DNA ( ccfDNA )  
Used for blood collection, anticoagulation, storage, transport, stabilization of ccfDNA

- 稳定期限**  
**Stability Period**  
 室温储存7天 ( 15-25°C )，超过35°C低于40°C不少于24小时。  
 7 days at room temperature (15-25°C), not less than 24 hours at 35-40°C
- 标识**  
**Identification**  
 浅紫色胶塞/透明安全头盖，标签印有不同编号的二维码  
 Light purple rubber stopper / transparent safety cap, unique number of 2-D barcode on label

产品代码 No.	材质 Material	采血管尺寸 Size	添加剂体积 Additive Volume	公称容量 Nominal Capacity
25022471	PET	13*100mm	0.5ml	4.5ml
20023471	PET	16*100mm	1ml	9ml



## ➤ Labtub® 血液cfRNA管 Labtub® Blood cfRNA Tube

用于采集、抗凝、储存、运输血液标本，稳定cfRNA

Used for blood collection, anticoagulation, storage, transport, stabilization of cfDNA

- 稳定期限

### Stability Period

室温储存7天 (15-25°C), 超过35°C不少于24小时

7 days at room temperature (15-25°C), not less than 24 hours at over 35°C

- 标识

### Identification

浅蓝色胶塞/透明安全头盖

Light blue rubber stopper / transparent safety cap

产品代码 No.	材质 Material	采血管尺寸 Size	添加剂体积 Additive Volume	公称容量 Nominal Capacity
25022491	PET	13*100mm	0.5ml	4.5ml
20023491	PET	16*100mm	1ml	9ml



## ➤ Labtub® 血液PRP管 Labtub® Blood PRP Tube

用于富含血小板血浆的采集分离。

Used for abbreviation of platelet rich plasma.



产品代码 No.	采血管尺寸 Size	血液抽吸体积 Blood Volume	PRP血浆提取量 PRP Plasma Volume
20013071	16×100mm	10ml	约为3.2ml
28033071	16×100mm	8ml	约为2.5ml
26033071	16×100mm	6ml	约为1.8ml
21219071	16×120mm	12ml	约为4ml
20039071	16×120mm	10ml	约为3.2ml
28039071	16×120mm	8ml	约为2.5ml
24039072	16×120mm (双层)	4ml	约为1ml
11314075	16×125mm (玻璃)	13ml	约为4.3ml
11134075	16×125mm (玻璃)	11ml	约为3.6ml
19034075	16×125mm (玻璃)	9ml	约为2.9ml
17534075	16×125mm (玻璃)	7.5ml	约为2.3ml
22036075	30×115mm	20ml	约为7ml
24036075	30×115mm	40ml	约为14.5ml
25036075	30×115mm	50ml	约为18ml

▶ **单个核细胞凝胶分离管(CPT管)**  
Single Nuclear Cell Gel Separation Tube (CPT Tube)

用于淋巴细胞免疫功能检测；HLA检测；残留白血病基因检测。  
一次离心直接分离淋巴细胞和单核细胞。

Cell immune function detection ; HLA detection ; Gene detection

产品代码 No.	材质 Material	采血管尺寸 Size	添加剂体积 Additive Volume	血液抽取体积 Blood Volume
14032065	玻璃(Glass)	13*100mm	0.1ml	4ml
18034065	玻璃(Glass)	16*125mm	0.1ml	8ml
18034065	玻璃(Glass)	16*125mm	135usp	8ml



▶ **ACD管 ACD Tube**

用于亲子鉴定、DNA检测及血液学分析ACD抗凝管。

Used for the paternity testing, DNA detection and hematology.



产品代码 No.	尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
15032151	13×100	PET	5ml	枸橼酸钠葡萄糖/ Sodium citrate glucose	1800支/箱
25032151	13×100	玻璃 (Glass)	5ml		1800支/箱

▶ **微量元素管 Trace Elements Tube**

专门用于微量元素的检测与分析。

Used for detection and analysis of trace elements.



产品代码 No.	尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
23051121	13×75	PET	3ml	肝素钠/Sodium Heparin	2400支/箱
25052121	13×100	PET	5ml		1800支/箱

### » PTH-rP专用采血管 (同位素专用采血管)

PTH-rP Tube (Isotope Dedicated Blood Collection Tube)

用于肿瘤标志物-甲状旁腺激素相关蛋白 (PTH-rP) 测定。

Used for tumor markers - Parathyroid Hormone Related Protein (PTH rP) detection .



尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
13×75	PET/玻璃(Glass)	2ml	EDTA-K2、抑肽酶 (Aprotinin)	2400支/箱
13×100	PET/玻璃(Glass)	2ml		1800支/箱

### » CTAD检测管 CTAD Tube

用于凝血因子测定，添加剂含有：枸橼酸钠；茶碱；腺苷和潘生丁，可以稳定凝血因子。

Used for detection of coagulation factor, Additive containing: citron acid sodium, theophylline, adenosine and dipyridamole, stabilize coagulation factor.



尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
13×75	PET	2ml	枸橼酸钠；茶碱；腺苷和潘生丁 Sodium citrate, theophylline, adenosine, dipyridamole	2400支/箱
13×100	PET	5ml		1800支/箱

### » RAAS专用采血管 RAAS Tube

用于肾素—血管紧张素—醛固酮 (RAAS) 测定 (高血压三项)。

Used for Renin - Angiotensin - Aldosterone (RAAS) detection (three hypertension).



尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
13×75	PET/玻璃 (Glass)	3ml	EDTA-K2,8-羟基喹啉二硫基丁二酸钠 (8-hydroxyquinoline, 2 thiol succinic, sodium)	2400支/箱
13×100	PET/玻璃 (Glass)	5ml		1800支/箱

### » 铅元素检测管 Lead Element Tube

专门用于进行铅元素检测与分析。

Used for detection and analysis of lead element.



尺寸 Size	材质 Material	规格 Volume	添加剂 Additive	每箱数量 Packaging
13×100	PET	5ml	肝素锂/ Lithium Heparin	1800支/箱