

静脉采血针 Blood Collection Needle

▶ 静脉采血针 Blood Collection Needle

静脉采血针与真空采血管配套使用，供医疗机构用于采集人体静脉血液样本。

Blood collection needle is matching use with blood collection tube, to collect human venous blood specimen for clinical routine detection by medical institutions.



项目	规格/Volume	每箱数量/Packaging
持针器/Holder		100个/盒
直针A/Pen Type A	21G	5000支/箱
直针B/Pen Type B	22G	5000支/箱
单翼A/Butterfly Type A	20~23G	2000支/箱
单翼B/Butterfly Type B	20~23G	4800支/箱



末梢采血产品

Capillary Blood Collection Product

▶▶ 末梢采血管 Micro Blood Collection Tube

用于采集末梢血液样本，添加剂种类有乙二胺四乙酸二钾（EDTA-K2）、肝素钠（肝素锂）、促凝剂等多种规格。

Collect capillary blood specimen, with additive of EDTA-K2, Sodium Heparin (Lithium Heparin), Clot Activator, etc.

添加剂	色标
乙二胺四乙酸二钾（EDTA-K2）	■
肝素钠（肝素锂）/Sodium/Lithium Heparin	■
促凝剂/Clot Activator	■



▶▶ 血液样本采集卡 Blood Collection Card

采用WHO推荐采血点应至少三个采血点的棉绒采样滤纸进行全血样本收集与保存，样本可保存18年。可用于多种生化、免疫学及DNA检测。

Collect and preserve the whole blood specimen by at least three blood collection points on the cotton velvet sampling filter paper, which is recommended by WHO, Specimen can be kept 18 years, used for a variety of biochemistry, immunology and DNA detection.



产品组成/Composition	特点/Feature
采集卡 /Acquisition Card	三个/ 四个/五个采血点/Three/Four/Five Blood Collection Points
酒精棉片/Alcohol Cotton	2片/2 Pills
采血针/Blood Collection Needle	21G /21G